





SeeBreath

CapnoDura Combi

CO₂ indicator/detector



Product Features

- * Clear indication of carbon dioxide helps verify endotracheal tube placement
- * Capable of indicating expired CO₂ levels continuously for up to 24 hours
- * Designed for adults and children
- * Quick response to exhaled CO₂ with a colour change from blue to yellow
- * Breath to breath visualisation up to 25 breath per minute

Specifications

Recommended for patients over 10 kg

Internal deadspace: 16 ml

Resistance to flow: 2 cm H₂O at 60 l/min flow

Connector ports: patient end : ISO 22mm OD/15 mm ID, circuit end : ISO 15mm OD

Shelflife: 2 years when stored at or below 30 degrees C