





SeeBreath

CapnoDura Pedi

CO₂ indicator/detector



Product Features

- * Clear indication of carbon dioxide helps verify endotracheal tube placement
- * Capable of indicating expired CO₂ levels continuously for up to 24 hours
- * Designed for pediatric use
- * Quick response to exhaled CO₂ with a colour change from blue to yellow
- * Breath to breath visualisation up to 40 breath per minute

Specifications

Recommended for patients between 1-15 kg

Internal deadspace: 3 ml

Resistance to flow: 2.0 cm H₂O at 10 l/min flow

Connector ports: patient end : ISO 15 mm ID, circuit end : ISO 15 mm OD

Shelflife : 2 years when stored at or below 30 degrees C